# **Getting Fewer Refractive Surprises.** ZEISS IOLMaster 700

2022-11-10 Camilla Arnesson Regional Sales Manager Nordics







• IOLMaster 500

• IOLMaster 700

# **Defining biometry for 20 years – ZEISS IOL Master** History of Optical Biometry



### **Passion & Commitment**

### Since 1999 ZEISS has

- Performed >800.000 simulated biometry test measurements
- Invested >1 million hours of development
- Released >50 software versions







### **Innovative Strengths**

### Key Innovations ZEISS IOL Master

- Telecentric Keratometry Automatic Measurement
- Composite Signals
- On-board post refractive surgery calculation
- First fully integrated markerless toric workflow

### Key Innovations ZEISS IOLMaster 700

- First SWEPT Source Biometry<sup>®</sup>
- Unique fixation check
- Patented Cornea-to-Retina Scan
- Up to 99% cataract penetration rate
- On-board toric IOL calculation
- Total Keratometry (TK<sup>®</sup>)
- Central Topography

## **Technology and Benefits** Speed and Precision



# This case study gives indication that **ZEISS** IOLMaster 700 cuts tech time by 60%.

#### Singh, CRSTE 2017

- 42 patients.
- ZEISS IOLMaster 500 vs, ZEISS IOLMaster 700
- Biometry from time the technician began entering data into the computer and instructing patients until the technician completed the scan on the second eye.





Check quality of measurements in the **'Quality check'** screen.

ZEISS

- Check quality of Central Topography maps by tapping on 'Topo'.
- Refer to the *IOLMaster* 700\_Plausability & Quality Check Guide\_EN\_32\_ 012\_00301 for further information on evaluation of measurements.

### **Technology and Benefits Overview**



# **ZEISS IOL**Master 700 with 3-zone telecentric keratometry

- Distance-independent Telecentric Keratometry
- Utilizes 18 telecentric points
- Shows excellent agreement with manual keratometry while achieving higher precision



Constant spot distance irrespective of device-to-eye distance

# **3-zone telecentric keratometry**





# **3-zone telecentric keratometry**





# Technology and Benefits Overview

**SWEPT Source Biometry Measurements** 



# **Biometry measurements are based on SWEPT SOURCE OCT images**

- Complex interpretation of graphs become obsolete
- Sources of errors caused by incorrect measurements and leading to refractive surprises can be eliminated



## **Technology and Benefits Overview** SWEPT Source Biometry Measurements





### **SWEPT Source Biometry**

Scan depth:	44 mm
Scan width:	6 mm
Fixation check width:	~1 mm (on retina
Resolution in tissue:	22 µm
Speed:	2,000 scans/s

#### **Standard Deviation of Repeatability**

Central cornel thickness	±2μm
Anterior chamber depth	±10 µm
Lens thickness	±19 µm
Axial length measurement	±9µm
Keratometry (SE)	±0.07 D

# **Technology and Benefits Overview** SWEPT Source Biometry





SWEPT Source Biometry<sup>®</sup> facilitates a Full-Length OCT Scan for visually verifying measurements and detecting unusual eye geometries



Fixation Check and Telecentric Keratometry for robust measurements

#### Patient Cölestin Cecile

Date of birth Patient ID	1/17/1954 000000014	Gender	Female	
Physician	Surgeon	Operator	Surgeon	

Date of calibration test:		by:		Result:	
Date of measurement:	3/1/2016	n:	1.3375	CVD:	12.00 mm

	right	Analyze				
			Eye status			
LS: Ref:	Phakic -2.25 D -0.75 D @ 34	VS: VA:	Vitreous body 20/25	1	LVC: LASIK	
			Biometric values	9		
AL: CCT. ACD: LT:	25.24 mm 489 μm 4.03 mm 3.59 mm	SD: 5 μm SD: 7 μm SD: 9 μm SD: 2 μm	WTW P	11.3 mm 5.2 mm	IX: +0.2 mm IV: - CW-Chord: 0.2 mm @	+0.0 mm 2 121°
SE: K1: K2: ΔK:	45.98 D (!) 44.84 D @ 173° 47.18 D @ 83° -2.34 D @ 173°	SD: 0.15 D SD: 0.75 D SD: 0.10 D	TSE: TK1: TK2: ΔTK:	37.76 D 37.66 D @ 37.87 D @ -0.21 D @	SD: 0.38 D 0° SD: 1.02 D 90° SD: 0.75 D 0°	
			B scan			
	and the second se					
	Keratometry		White-to-	-white	Fixati	ion
		r		0		
		Velue une offici	manually No un	haruseem aul	M NEWGA	



# Technology Deep Dive

**Dense Cataracts Penetration** 



# Rate of successful scans with ZEISS IOLMaster 700 was estimated to be 99.5 %

#### Hirnschall et al. Ophthalmol Ther, 2018

- 1,226 eyes of 613 patients (23 unsuccessfully measured eyes included into the study)
- Biometrical examination with ZEISS IOLMaster 500 and ZEISS IOLMaster 700
- 91.3 % (21/23) of the eyes that were measured unsuccessfully with ZEISS IOLMaster 500 were measurable with the ZEISS IOLMaster 700
- The estimated overall rate of unsuccessful scans with the SS-OCT device was 0.5% (6/1226)



Successful (top) and unsuccessful (bottom) longitudinal SS-OCT scan. The arrow depicts the macula

**Hirnschall N, Varsits R, Doeller B, et al.** Enhanced Penetration for Axial Length Measurement of Eyes with Dense Cataracts Using Swept Source Optical Coherence Tomography: A Consecutive Observational Study. Ophthalmol Ther 2018;126(4):524

### **Technology Deep Dive** Dense Cataracts Penetration

ZEISS

# Cataract penetration rate of up to 99%, Reduction in need for ultrasound of up to 92%

#### Varsits et al, ESCRS 2016

- The ZEISS IOLMaster 700 with SWEPT Source OCT is designed to optimize workflow efficiency, even when handling dense cataracts.
- A comparative clinical study with more than 1,200 eyes



**Image: Prof de la Torre Estremadoyro M.** Case Study: IOLMaster 700 in extreme cataract. CRSTE 2015:700

# Technology Deep Dive

**Dense Cataracts Penetration** 



# In this case study ZEISS IOLMaster 700 could measure axial length in very dense nucleous cataracts

#### De La Torre, CRSTE 2015

- 2 Patients with very dense nucleous cataracts, including brunescent, rubra, nigra and white
- axial length measurements with previous optical biometry devices weren't possible



Fig. 1: IOLMaster 700 OCT image with measurement calipers Fig. 2: White instumescent cataract Fig. 3: Slit lamp image of the cataracta ruba eye

Prof de la Torre Estremadoyro M. CASE STUDY : IOLMASTER 700 IN EXTREME CATARACT. 2015:700.





# Fixation check will allow you to detect poor fixation during measurement

- 1 mm Foveal scan helps verify if the patient has fixated correctly
- Incorrect fixation can lead to incorrect axial length measurement and wrong keratometry data and thus refractive surprises

Helps to reduce the risk of refractive surprises



# Example: More accurate & repeatable measurements





# Example: More accurate & repeatable measurements



# Example: More accurate & repeatable measurements





# Example: helps to identify measurement errors





# Example: helps to identify measurement errors



# Example: helps to identify measurement errors



ZEINS

### **Example: Influence on K-values**





1D difference in calculated IOL power!



# In this case of Srivannaboon a macular hole could be observed with ZEISS IOLMaster 700 before cataract surgery was performed

#### Srivannaboon, CRSTE 2015

- SS-OCT enabled observation of abnormal macula morphology\* in this case
- Additional OCT examination was needed



Fixation check of ZEISS IOLMaster 700 overlaid CIRRUS 4000 retina OCT

Srivannaboon, Case Study: IOLMaster 700, CRSTE, June 2015

\* Findings need to be verified and pathologies diagnosed with a dedicated retina OCT or other clinical standard methods.



**Regular Fixation Check** 







Figure 1. SWEPT Source Biometry of our patient with the IOLMaster 700 showing intraretinal fluid (Fixation Check image on the right)\*

- The unique Fixation Check supports to detect poor patient fixation
- It may also help to indicate unusual eye structures for better patient selection<sup>1</sup>
- An incidential finding (e.g. BRVO as shown here) may prompt for a comprehensive OCT examination.



Figure 2. SD-OCT image of the left eye of our patient\*

1 As the ZEISS IOLMaster 700 is clearly not intended to be used for diagnostics, findings need to be verified and pathologies diagnosed with a dedicated retina OCT ate, Image courtesy of Prof. O. Findl, Hanusch Hospital Vienna, Austria



Retinal Pathology: Macular Foramen



Image courtesy of Prof. W. Sekundo, Philipps University Hospital Marburg, Germany



Retinal Pathology: Macular Pucker



Image courtesy of Prof. W. Sekundo, Philipps University Hospital Marburg, Germany



Retinal Pathology: AMD, RPE detachment



Image courtesy of Prof. W. Sekundo, Philipps University Hospital Marburg, Germany

# **Technology Deep Dive**

**Total Keratometry** 





## **Total Keratometry (TK)** Replacing assumptions with measurements



- **Total Keratometry(TK**®) is a new measurement that combines telecentric keratometry and SWEPT Source OCT technology for the assessment of anterior and posterior corneal curvature.
- The purpose of TK is to replace standard keratometry, aiming to help to reduce outliers and improve refractive outcomes of IOL calculation in cataract surgery.
- TK is ULIB-compatible, therefore existing standard formulas and IOL constants may be applied.



Picture source: Carl Zeiss Meditec media database

## **Remarkable IOL Calculation Results Observed**

Improve toric IOL calculation



Figure 3: Outcomes of toric IOL calculations with the Haigis-T formula. CYL APE: Absolute prediction error for cylinder; frequency of eyes in respective CYL APE diopter ranges; N=145 eyes\*. Figure 4: Outcomes of toric IOL calculations with classic Barrett Toric Calculator and the new Barrett TK Toric formula; CYL APE: Absolute prediction error for cylinder; frequency of eyes in respective CYL APE diopter ranges; N=145 eyes\*.

#### Source: Total Keratometry Compendium

Fabian and Wehner 2018

\* \*Retrospective post-hoc analysis of 145 normal cataract eyes implanted with aspheric IOL, 6 weeks post-op.

# Improve toric IOL calculation.



# **Total Keratometry (TK)**

Replacing assumptions with measurements



*"I don't recall ever seeing such a dramatic improvement, by simply adding an additional parameter."* 

"Despite the minimal change (in the formula) the spherical prediction is definitely improved with BUII-TK"



**Graham Barrett** 

# In this study in post

In this study in post-myopic LASIK eyes, Barrett True-K with TK improved the outcome prediction compared to Barrett True-K with classic K's within ±0.5 D by >12% (p = 0.04)

#### Lawless et al., Clin Experiment Ophthalmol 2020

- 72 eyes of 50 patients
- The Barrett True-K TK provided the lowest mean refractive prediction error and variance for both prior myopes and hyperopes undergoing cataract surgery





## Total Keratometry (TK)

Replacing assumptions with measurements

# Technology Deep Dive Central Topography

- Toric and multifocal IOL implantation requires exclusion of irregular corneas.
- Corneal asymmetries in a central optical zone of 4-5 mm are regarded as clinically relevant in this regard.
- With no changes in workflow the ZEISS IOLMaster 700 measures central corneal topography.
- It provides anterior and total axial power maps, designed to detect visually relevant central corneal asymmetries.

### Total Power Map

Pupil Diame Topography Step	ter Diameter	3 mm 4.8 mm 0.5 dpt	
49.5	44.0	38.5	





# Technology Deep Dive Central Topography





## Connectivity Deep Dive ZEISS EQ WORKPLACE



# With ZEISS EQ Workplace you get access to your biometry data anytime, anywhere

The new "heart" of the ZEISS Cataract Workflow

- Remotely calculate IOLs and plan surgeries, facilitating the exchange of data between diagnostics and the OR.
- Connect the ZEISS IOLMaster and other devices in the cataract workflow to ZEISS FORUM and ZEISS CALLISTO eye, allowing a seamless toric and premium cataract workflow.
- A new level of protection against never-events: ZEISS IOLMaster data is transferred automatically to the ZEISS EQ Workplace, to populate the relevant fields for IOL calculation and selection.



# Connectivity Deep Dive ZEISS CALLISTO eye

# Save time, increase efficiency and reduce residual astigmatism

#### Computer assisted cataract surgery

- Manual marking steps can be skipped altogether for an efficient and precise toric IOL alignment to reduce residual astigmatism.
- Starting with a biometry reference image from the ZEISS IOLMaster, data is transferred smoothly to ZEISS CALLISTO eye. This data is used to create overlays in the eyepiece.





# Seeing beyond